

# Cherishmet Inc.

301 Route 17, Suite 800, Rutherford, NJ 07070 USA  
Web site: <http://www.cherishmet.com>



Telephone: (201) 842-7612  
Facsimile: (201) 455-6816  
Toll Free: (866) 804-0083

## AG- 900, AG-1000

Model: AG-900, AG-1000

### Application:

AG-900, AG-1000 activated carbon is a virgin activated carbon. Produced from bituminous coal. It is granular in form. Its ideal for many liquid phase applications including the removal of organics from water stream. Its superior adsorptive capacity and surface area make it particularly attractive for use with potable water while its size make it appealing where low-pressure drop is need.

### Specification:

	ASTM TEST METHODS
Mesh Size:	6x10, 6x12, 8x16, 8x20, 8x30, 12x40, 4x10
Ash(WT%):	18 max. D2866-94
Moisture:	3% max. D2867-91
Hardness No.	90 min. D3802-79
Iodine No.	900mg/g min. D4607-94
Density:	0.45-0.52g/ml D2854-94
Abrasion:	80min. B604

### Packaging:

Material to be packed in black super sacks with 500kg or 1100lb net each, 2 sacks per pallet. 20-25kg small bags are available.

### Remarks:

The above specification can be adjusted in accordance with the customer's requirement.